


# BARRIER SYSTEM FOR SURGICAL THEATRES

## Sterile Surgical Drapes and Sets

### Dentistry


#### Dental Set 801-001-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	60 x 60	1
Instruments table drape	2SBL	75 x 90	1
Dental drape with aperture	DS3	150 x 150	1
Surgical drape	DS3	45 x 75	2
OP tape	AT+Vis	10 x 50	2
Cord cover	PE	8 x 120	2



#### Surgery Set - Implantology 801-013-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 80	1
Instruments table drape	2SBL	80 x 120	1
OP tape	AT+Vis	5 x 20	2
Surgical drape	SMS	75 x 90	1
Dental drape with adhesive aperture	SMS	120 x 265	1
MAYO cover	PE/DS2	80 x 145 / 50 x 145	1
Absorbent towel	Airlaid cellulose	30 x 33	2
Flat bowl	PP	125 ml	1
Non-woven compress, 4ply	NW	5 x 5	3
Twisted tampon	Gauze	15 x 15	3
Surgical gown Prima	SMS	L	2



#### Dental Set 60-19-126-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	100 x 100	1
Dental drape	DS3	150 x 200	1
Tray 4 compartment	PP	15 x 22	1
Yankauer suction handle	PP		1
Gauze compress	Gauze	10 x 10	10
Cord cover	PE	8 x 120	2
Light handle cover with adhesive edge	SMS	15 x 20	4
Surgical drape with adhesive edge	DS2	75 x 75	3
Absorbent towel	Airlaid cellulose	30 x 33	2
Wrapping sheet	FS2	75 x 80	1
Surgical gown Prima	SMS	L	2
Cap Baret Standard	SB		2
Cap Forest Standard	SB		1

