


BARRIER SYSTEM FOR SURGICAL THEATRES

Sterile Surgical Drapes and Sets

Gynaecology


Gynaecology Set 601-103-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 80	1
Instruments table drape	2SBL	100 x 150	1
Reinforced drape with adhesive edge	2SBL+Vis28	180 x 200	1
Perineum drape, adhesive aperture	DS2	75 x 120	1
Limb cover	SMS	80 x 120	2
Surgical drape	2SBL	100 x 150	1
Surgical drape	SMS	50 x 50	1
OP tape	AT+Vis28	10 x 50	2
MAYO cover	PE/SB28	80 x 145/50 x 145	1
Surgical gown Classic	SMMMS	size L	1
Absorbent towel	Airlaid cellulose	30 x 33	4




Gyn Set for Laparoscopy 601-132-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 80	1
Instruments table drape	2SBL	150 x 200	1
Absorbent towel	Airlaid cellulose	30 x 33	2
Surgical drape	DS2	75 x 90	1
Limb cover	SMS	80 x 120	2
Surgical drape with adhesive aperture and edge	2SBL+Vis28	180 x 200	1
MAYO cover	PE/SB28	80 x 145/50 x 145	1



LAVH Set 601-107-S


Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	100 x 100	1
Instruments table drape	2SBL	100 x 150	1
Perineum drape, adhesive aperture, pouch	DS2	75 x 120	1
Limb cover	SMS	80 x 140	2
Laparotomy drape reinforced, IF	2SBL+Vis28	180 x 200	1
Surgical drape	SMS	75 x 120	3
MAYO cover	PE/SB28	80 x 145/50 x 145	2
Surgical drape	2SBL	100 x 150	1
Camera cover	PE	14 x 250	1
Syringe LS	PP	10ml	1
Scalpel blade type 20			1
Scalpel blade type 11			1
Tampon twisted, RTG	Gauze	30 x 30	40
Tampon twisted, RTG	Gauze	20 x 20	20
Gauze compress	Gauze	10 x 20	4
Folded gauze	Gauze	15 x 30	4
Surgical gown Classic	SMMMS	XL	3
Surgical gown Classic	SMMMS	L	1



Gynaecology


Caesarean Set 602-153-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	100 x 100	1
Instruments table drape	2SBL	100 x 150	1
Anaesthetic drape	SMS	100 x 155	1
Caesarean drape with PUR brace pouch and IF	DS2	200 x 210	1
OP tape	AT+Vis28	10 x 50	1
MAYO cover	PE/SB28	80 x 145 / 50 x 145	2
Surgical gown Classic reinforced	SMMMS	XL	1
Surgical gown Classic	SMMMS	XL	2
Surgical gown Classic	SMMMS	M	1



Caesarean Set 602-121-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 90	1
Caesarean drape with IF	DS2	240 x 320	1
MAYO cover	PE/SB28	80 x 145 / 50 x 145	1



Caesarean Set 602-147-S

Component	Material	Size (cm)	Pcs
Instruments table drape	FS2	150 x 190	1
Caesarean T-shape drape, PUR brace pouch and IF	2SBL	180 x 320	1
Absorbent pad Santex	Cellulose	60 x 90	1
Absorbent towel	Airlaid cellulose	30 x 33	2
MAYO cover	PE/SB28	80 x 145 / 50 x 145	1

