



# Ophthalmology

## Eyelid Pack

**40-14-028-S**

Component	Material	Size (cm)	Pcs
Instruments table drape	2SBL	100 x 150	1
Ophthalmic drape with adhesive aperture, 2 pouches	SMS	95 x 100	1
Transparent bowl	PE	250 ml	1
Twisted tampon	Gauze	30 x 30	4
Hypodermic needle grey		0,4 x 40, 27G x 1 1/2"	2
Hypodermic needle pink		1,2 x 40, 18G x 1 1/2"	1
Non-woven compress, 4ply	NW	7,5 x 7,5	10
Syringe LL	PP	2ml	2




# BARRIER SYSTEM FOR SURGICAL THEATRES

## Sterile Surgical Drapes and Sets

### Ophthalmology


#### Set PHACO 201-083-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 80	1
Ophthalmic drape with IF, 2 pouches	SMS	120 x 160	1
Surgical drape with adhesive edge	SMS	52 x 80	1
Handle cover	SMS	20 x 50	2
Device cover PHACO	SMS	42 x 57	1
Instruments table cover, reinforced	SMS, SBPE	50 x 77,5	1



#### Cataract Set 201-132-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	DS2	50 x 70	1
Handle cover	SMS	20 x 75	2
IF	PU	15 x 20	1
Surgical drape with adhesive aperture	DS2	150 x 180	1
Examination gloves	Latex	size 7-8	1
MAYO cover	PE/SB28	80 x 145/50 x 145	1
Drying tampon	NW+ cellulose	0,5 x 6,5	5
Hypodermic needle brown		26G x 1/2", 0,45 x 12 mm	1
Syringe LL	PP	10ml	1
Syringe LL	PP	5 ml	1
Syringe LS	PP	2 ml	1
Insulin syringe U-100, with needle	PP	1 ml / 30G	1
Twisted tampon	Gauze	30 x 30	2



#### Ophthalmic Set 60-04-034-S

Component	Material	Size (cm)	Pcs
Instruments table drape	FS2	60 x 60	1
Ophthalmic drape with perforated IF	SMS	50 x 60	1
Non-woven compress, 4ply	NW	10 x 10	5
Gallipot bowl	PP	60 ml	2
Cotton swab, small bud		15	2
Double ended scleral marker, flat	PP	0,35 - 0,40	1
Eye speculum (solid blades)	Stainless steel		1
Kocher	Stainless steel	13	1

