


BARRIER SYSTEM FOR SURGICAL THEATRES

Sterile Surgical Drapes and Sets

General Surgery


Laparotomy drape with AT edge and IF.

Code	101-116-S
Material	DS3
Size (cm)	240 x 320
Reinforcement (cm)	
Aperture (cm)	25 x 30
Accessories	IF, AT




Laparotomy drape, reinforced, with AT edge and IF, short-term surgeries.

Code	101-110-S
Material	2SBL
Size (cm)	200 x 320
Reinforcement (cm)	Vis28 100 x 150
Aperture (cm)	25 x 30
Accessories	IF, AT




Laparoscopy drape, reinforced, with IF and leggings.

Code	105-130-S
Material	SMS
Size (cm)	260 x 280
Reinforcement (cm)	DS2 80 x 150
Aperture (cm)	25 x 30
Accessories	IF, leggings



Laparotomy drape with IF, short-term surgeries.

Code	101-207-S
Material	FS2
Size (cm)	200 x 320
Reinforcement (cm)	
Aperture (cm)	25 x 30
Accessories	IF



Laparotomy drape with adhesive aperture, short-term surgeries.

Code	101-208-S
Material	FS2
Size (cm)	200 x 320
Reinforcement (cm)	
Aperture (cm)	25 x 30
Accessories	AT

