


# Urology


Combined TUR drape, adhesive aperture with collection pouch.

Code	603-205-S
Material	DS3+PE
Size (cm)	220 x 300
Aperture (cm)	6 x 10
Accessories (cm)	AT, pouch




Combined TUR drape with aperture.

Code	603-126-S
Material	SMS+DS2
Size (cm)	(150 x 160) + (50 x 110)
Aperture (cm)	5,5 x 12,5



Urology drape reinforced, with 2 apertures and collection pouch.

Code	603-209-S
Material	SMS
Size (cm)	200 x 260
Reinforcement (cm)	DS2 30 x 70
Aperture (cm)	5 and 6 x 8
Accessories (cm)	Pouch



Surgical drape with adhesive edge and adhesive aperture.

Code	604-108-S
Material	2SBL
Size (cm)	180 x 200
Reinforcement (cm)	Vis 75 x 100
Aperture (cm)	20 x 2 x (22,5)
Accessories (cm)	AT

