


# BARRIER SYSTEM FOR SURGICAL THEATRES

## Sterile Surgical Drapes and Sets

ORL


### ORL Set 106-056-S

Component	Material	Size (cm)	Pcs
Instruments table drape	FS2	150 x 160	1
Surgical drape with adhesive aperture	FS2	170 x 240	1
Absorbent towel	Airlaid cellulose	30 x 33	1



### Head and Neck Surgery Set 106-064-S

Component	Material	Size (cm)	Pcs
Instruments table drape	FS2	150 x 190	1
U-drape, adhesive opening 6 x 40 cm	FS3	150 x 125	1
Surgical drape with adhesive edge	FS3	200 x 200	1
Absorbent towel	Airlaid cellulose	30 x 33	4
MAYO cover	PE/SB28	80 x 145 / 50 x 145	1



SETS  
SETS

