BARRIER SYSTEM FOR SURGICAL THEATRES

Set's Accessories

Surgical gowns

Properties:

Surgical gowns provide high protection from body fluids as well as infection.

Gowns are made from non-linting nonwoven hydrophobic material SMMMS or SPUNLACE, with antistatic treatment. Both materials are breathable and comfortable to wear.

Gowns are of a tie-round type and fastened by the four tie system to keep the wearer aseptic.

Sewn-in sleeves are finished with knitted elastic cuffs, the upper seam of sleeves is done by ultrasonic welding, as well as ties.

Surgical gowns with reinforcement – are equipped with a liquid proof material SBPE at both the front side and on sleeves. Although the material is laminated, it is breathable and provides high level of protection.

Neck Colour Differenciation indicates Gowns' liquid resistance:

 Blue
 low fluid contact surgery

 Yellow
 middle fluid contact surgery

 Green
 middle and high fluid contact surgery

 Red
 long-lasting and high fluid contact surgery



PRIMA


Made of hydrophobic non-woven textile SMS.

Size	PRIMA Ref	PRIMA Reinforced Ref	Length (cm)
M	354523-S	354623-S	115
L	354533-S	354633-S	125
XL	354543-S	354643-S	135
XXL	354553-S	354653-S	145
XXL-XL	354563-S	354663-S	160

CLASSIC

Made of hydrophobic non-woven textile SMMMS.

Size	CLASSIC Ref	CLASSIC Reinforced Ref	Length (cm)
M	351023-S	351123-S	115
L	351033-S	351133-S	125
XL	351043-S	351143-S	135
XXL	351053-S	351153-S	145
XXL-XL	351063-S	351163-S	160



PREMIUM

Made of hydrophobic non-woven textile SPUNLACE. This type of gowns is designed for prolonged surgeries. Very comfortable to wear.

Size	PREMIUM Ref	PREMIUM Reinforced Ref	Length (cm)
M	350023-S	350123-S	125
L	350033-S	350133-S	135
XL	350043-S	350143-S	145
XXL	350053-S	350153-S	155