

DISPOSABLE CLOTHING FOR HOSPITALS

Surgical Face Masks

- Surgical Face Masks Dina Standard Type II-white, with active silver particles in the filtration layer, non-sterile
- Bacterial Filtration Efficiency >99,5 %
- Viral Filtration Efficiency >99,6 %
- Tested in Nelson Laboratories, USA
- Made in Czech Republic
- Latex free
- Conform to EN 14683 Typ II standards
- CE

Surgical Face Masks Dina-Hitex are primarily used in hospitals and other medical facilities. At present, they are also in great demand for common use as a prevention against droplet transmission of infections.

Thanks to the highly efficient filtration layer, DINA-HITEX Face Masks provide maximum protection (more than 99%) against the penetration of bacteria and viruses and at the same time provide high comfort even when worn for longer periods. FM do not irritate the skin and are pleasant to breathe.

The integrated shaping wire adapts well to your face shape.

DINA-HITEX Face Masks are intended for single use.

Place a clean FM on the face with a shaping wire upwards and shape it around the root of the nose with gentle pressure. Attach the ties around the head or the rubber band behind the ears and then spread the folds so that they completely cover the area of the nose and mouth.

Dispose of the used Face Mask.

Standard

Ref	Features	Packaging
20251 DH	Three-layer surgical FM with ties Thanks to the fixation around the head, it holds well during long-term surgical procedures.	50pcs/box, 500pcs/carton
20254 DH	Three-layer surgical FM with elastics The FM is fixed behind the ears and, thanks to its ease of use, is popular not only in hospitals, but also in gastronomy, cosmetics and common use during a day.	100pcs/box, 1000pcs/carton



Standard-Antireflex

Ref	Features	Packaging
20265 DH	Three-layer surgical FM ANTIREFLEX with ties Relieves the burden of vision in strong lighting in operating rooms. Recommended when using glasses. Thanks to the fixation around the head, it holds well during long-term surgical procedures.	50pcs/box, 500pcs/carton



Standard-Shield

Ref	Features	Packaging
23254 DH	Three-layer surgical FM with elastics and a shield Provides eye and upper face protection.	25pcs/box, 100pcs/carton



Goggles

Ref
92002

