


Liquid Collection Pouches

Adhesive Pouches

Designed to collect fluid from the surgery area.
Pouches are made of strong PE transparent foil.
Equipped with adhesive tape at the back side which ensures easy fastening to the drape.
Front part of the pouch is equipped with a forming brace.


1-Compartment Pouch

Size (cm)	Ref
22 x 29	60-0101-S
28 x 44	60-0102-S
28 x 55	60-0103-S




2-Compartment Pouch

Size (cm)	Ref
25 x 28	60-0124-S
28 x 44	60-0125-S



3-Compartment Pouch


Size (cm)	Ref
28 x 44	60-0149-S



Pouch for Neurosurgery

Adhesive tape and semi oval aperture on upper rear side edge enables perfect adjustment to the head shape.
The front part with forming brace enables free fluid access into the pouch.


Size (cm)	Ref
30 x 32	60-0111-S



Pouches with Filter

The rear side of the pouches is equipped with a 5 cm wide adhesive tape.
The front part of the pouch is equipped with a forming brace.
Conical drainage port for tube setting (6 - 10 mm in diameter) is equipped with a filter to prevent from clogging.

Ref	Size (cm)	Accessories
65-0101-S	50 x 50	Drainage, filter
65-0102-S	50 x 60	Drainage, filter
65-0106-S	60 x 90	Drainage, filter



BARRIER SYSTEM FOR SURGICAL THEATRES

Set's Accessories

Liquid Collecting Pouches

Arthroscopic Pouches - Knee, Shoulder

Pouches are designed to collect fluid during knee or shoulder arthroscopy. They are made of strong PE foil, equipped with a conical drainage port for tube setting (6 - 10 mm in diameter) and a forming brace. Tight shape around the limb is ensured by 2 flexible elastic membranes with oval aperture 6 x 8 cm.

Ref	Size (cm)	Accessories
60-0122-S	60 x 100	2 elastic membranes



Caesarean Pouches

Made of strong transparent PE foil. Pouches are equipped with an adhesive tape around the aperture. The front part has either a forming brace or a PUR foam brace to retain fluids safely. Pouches are accessoried with a conical drainage port for tube setting (6 - 10 mm in diam.) on both sides or with a sleeve outlet.

Ref	Size (cm)	Accessories
60-0106-S	80 x 100	PUR foam
60-0107-S	80 x 100	Plastic forming brace
60-0145-S	110 x 110	PE sleeve outlet



Pouches for Gynaecology

They are made of strong PE foil as either 1-compartment or 2-compartment pouches. All pouches are equipped with a forming brace in front and an adhesive tape on rear side to be attached to the drape. Pouches are fitted with a conical drainage port or sleeve outlet and volume indicator.

Ref	Size (cm)	Accessories
60-0146-S	60 x 150	1 compartment, without scale, sleeve outlet
60-0153-S	60 x 100	1 compartment, with scale and drainage port
60-0148-S	60 x 100	2 compartments, with scale and drainage port